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**The Automatic Ash-Colour-Turbidity Analyser measures:**

- The conductivity in accordance with the official ICUMSA method
- The colour of the filtrated sample solution
- The unfiltrated colour of the same solution in order to calculate the turbidity as the difference of the two absorption measurements

### Preparation time 1 minute!

This system saves your valuable time and makes your quality control much more efficient and reliable.

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### Technical Data

Colour in solution accuracy:  $\pm 1$  IU (up to 30 IU)  
Turbidity accuracy:  $\pm 10$  %  
Ash content accuracy:  $\pm 0.001$  %  
Accuracy Brix solution:  $\pm 0.1$  Bx

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## Automatic Ash-Colour-Turbidity Analyser

Inorganic components in sugar remaining after complete incineration of organic substances are defined as ash. As a purity criterion for sugar and other products of the sugar industry the conductometric determination of ash is preferred, because it is closely related to sulphated ash. Conductivity can be rapidly measured with reproducible results and is contained in most national and international sugar quality regulations.

The **Automatic Ash-Colour-Turbidity Analyser** of SCHMIDT+HAENSCH measures the conductivity in accordance with the official ICUMSA method in a 28 Bx solution (ICUMSA 2002, GS2/3-17). Also the colour is an important criterion for the quality of refined sugar products. According to ICUMSA method described in GS 2/3-10, 2002, the solution colour is measured in a filtrated 50 Bx solution in a 100 mm cuvette photometrically at 420 nm and given in ICUMSA colour units (IU). The system also measures the unfiltered colour of the same solution in order to calculate the turbidity as the difference of the two absorption measurements. This value is especially important in soft drink production as higher values will cause non-pleasant optics of the final product.

By combining these three important determinations in an automatic laboratory system SCHMIDT+HAENSCH helps to improve the daily work in sugar laboratory. The obtained results are automatically transmitted using laboratory data acquisition software realising a complete solution from one hand.

## System Description

The analyser is a stand alone PC-controlled system containing a colorimeter, a scale and a conductometer with all necessary pumps and valves for fluid handling needing only a few preparations before automatically analysing the sugar product:

- Fill in a certain amount of sugar
  - Insert a new filter
  - Close the door and start
  - All of the above will take only **one minute** to finish
- After ca. 25-30 minutes (depending on sample) a full analysis will be finished and a new run can be started immediately. System connection to be provided by the user:
- Distilled water
  - Compressed air (4 to 8 bar)

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